



ONCOLOGY

In vitro diagnostic tests for oncologic markers detection.

EXPERTEAM SRL

Experteam is an Italian company with an history of more than 20 years in developing, producing and distributing *in vitro* diagnostic products

OUR OFFER

Diagnostic solutions for Azoospermia, Fragile X Syndrome, Huntington Disease, Hemochromatosis, Microsatellite Instability, B and T-cell Lymphomas, HPV and Food Intolerance (Celiac Disease and Lactose Intolerance) among others.

RQ Assays, based on Multiplex real-time PCR or Allelic Discrimination Assays, are used to detect specific alleles or Single Nucleotide Polymorphisms (SNPs).

FL Assays, based on PCR amplification followed by Automated Capillary Electrophoresis.

MICROSATELLITE INSTABILITY

Code	Kit name	Technical specifications	# Tests	Label	Info
MS.01FL	HNPCC Kit 1 - FL	Capillary Electrophoresis Kit	40	CE/IVD	(for Applied Biosystems and Beckman Coult. genetic analysers)

EXPERTEAM srl

Via della Libertà 12 - 30175 Marghera Venezia

Tel. +39 041 509 3101 -Fax +39 041 509 3102

www.experteam.it

Commercial information Italy:

commerciale@experteam.it

Foreign trade information:

sales@experteam.it

Italian and foreign technical support: supportotecnico@experteam.it

www.experteam.it

CLONALITY AND LYMPHOMAS

Code	Kit name	Technical specifications	# Tests	Label	Info
BL.01FL	B cell lymphoma kit-FL	Capillary Electrophoresis Kit	40	CE/IVD	(for Applied Biosystems genetic analysers)
TL.01FL	T cell lymphoma kit-FL	Capillary Electrophoresis Kit	40	CE/IVD	(for Applied Biosystems genetic analysers)

CALRETICULIN

Code	Kit name	Technical specifications	# Tests	Label	Info
MP.01FL	MPNs kit 1-FL	Capillary Electrophoresis Kit	40	CE/IVD	(for Applied Biosystems genetic analysers)

PAPILLOMA VIRUS

Code	Kit name	Technical specifications	# Tests	Label	Info
PL.04RQ	HPV HR/LR 23types detection kit-RQ	Real Time PCR Kit	25	CE/IVD	HPV16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68, 6, 11, 26, 53, 67, 70, 73, 81,82

THIOPURINES

Code	Kit name	Technical specifications	# Tests	Label	Info
TP.01FL	TPMT variants kit-FL	Capillary Electrophoresis Kit	40	CE/IVD	c.238G>C; c.460G>A; c.719A>G (for Applied Biosystems genetic analysers)

CHROMOSOMAL TRANSLOCATION

Code	Kit name	Technical specifications	# Tests	Label	Info
To922SRQ	Translocation kit t(9;22) screening - RQ	Real Time PCR Kit	40	CE/IVD	BCR-ABL M-bcr, m-bcr, u-bcr screening