



PHARMACOGENETICS

IL28B

It has been recently demonstrated by studies of Genome Wide Association (GWA) that the single nucleotide polymorphisms (SNPs) rs12979860 and rs8099917 in the gene interleukin-28B (IL28B) coding for IFN- λ 3 are closely associated with virologic response for the treatment PegIFN-alpha-2a / RBV.

Experteam offers a complete system for allelic discrimination (genotyping) of rs12979860 and rs8099917 in the interleukin-28B (IL28B) gene.

Code	Name of kit	Technical specs	Amount of tests	Info	Additional info
CH.01RQ	CHC 28B Kit 1 - RQ	Real Time PCR Kit	40	CE/IVD	Allelic discrimination



CHC 28B Kit 1 – RQ: schematic results after allelic discriminations

