


## HPV (Real Time PCR)

Complete systems for screening of Papilloma virus and evaluation of high risk group (HPV 16, 18, 31, 33, 35, 52, 58) by Real Time PCR (cod. PL.02RQ) and for high risk genotyping of HPV 16,18,31,33,35,39,45, 52,58,59,66,68 (cod. PL.03RQ).

Code	Name of kit	Technical specs	Amount of tests	Info	Additional info
PL.02RQ	Papilloma Virus Kit 2-RQ	Real Time PCR Kit	50	CE/IVD	melting analysis
PL.03RQ	Papilloma Virus Kit 3-RQ	Real Time PCR Kit	50	CE/IVD	genotyping
mlg.hpv.50	F-HPV Typing Multiplex Fluorescent Kit	Capillary Electrophoresis Kit	50	CE/IVD	genotyping



Dissociation curve shows the internal control ( $\beta$  Globin) at 80° C and samples positive for HPV at 77° C

## Dot Blot Strip

Code	Name of kit	Technical specs	Amount of tests	Info	Additional info
VR-124.20.40	AlphaStrip II HPV	Dot Blot Strip	20	CE/IVD	Genotyping
VR-125.20.40	AlphaStrip III HPV	Dot Blot Strip	20	CE/IVD	Genotyping

